

ANKLIZER II

FEATURES

- Deluxe foam liner designed for greater comfort.
- Designed for stable fractures and severe ankle sprains.
- Rocker sole and low profile heel facilitates a more natural gait pattern.
- Wider foot bed provides greater comfort.
- Uprights are contoured for a more anatomical fit at the calf.
- Cushioned heel for greater energy absorption.



ITEM #	FEMALE SHOE SIZE	MALE SHOE SIZE
BC0279-XS	4 1/2	3 1/3
BC0282-S	5 TO 8 1/2	4 TO 7
BC0284-M	9 TO 11	7 1/2 TO 10
BC0285-L	11 1/2 TO 13	10 1/2 TO 12
BC0286-XL	13 1/2 +	12 1/2 +