

SWEDE-O ATOM

FEATURES

- Hinged ankle brace fits in most athletic shoes.
- Easy application and single strap adjustment.
- Patented footplate cleats eliminate shifting of brace during use.
- Airflow EVA side pads keep ankles cooler and drier.



ITEM # BLACK
ATOM-S
ATOM-M
ATOM-L

MENS SHOE
SMALL (5 - 8 1/2)
MEDIUM (9 - 12)
LARGE (12 1/2 - 16)

FEMALE SHOE
SMALL (6 1/2 - 9 1/2)
MEDIUM (10 - 13)
LARGE (13 1/2 - 16)