

BAKER - NIGHT SPLINT

FEATURES

- Designed to provide relief from pain and discomfort of plantar fasciitis and Achilles tendonitis.
- Allows the patient to apply consistent stretch to the plantar fascia by holding foot in dorsiflexion.
- Soft padded cover provides greater patient comfort.
- Removable foam wedge insert included for additional 5 degree stretch.
- Fits either left or right.



ITEM #	WOMENS SHOE SIZE	MENS SHOE SIZE
BC4802-S	5 1/2 TO 6 1/2	5 TO 6
BC4802-M	7 TO 8	6 1/2 TO 7 1/2
BC4802-L	8 1/2 - 9 1/2	8 TO 9
BC4802-XL	10 TO 11	9 1/2 TO 10 1/2
BC4802-XXL	11+	11+