

# BODYFLEX II KNEE BRACE

## FEATURES

- Universal J Buttress can be used as patella subluxation brace for a right or left.
- Can also be used for Chondromalacia and Jumpers knee.
- Movable straps for optimum placement.
- Seamless Body Flex back panel and quadriceps.
- 4 Spiral stays for medial and lateral stability.
- Patented silicone non-slip suspension band keeps brace in place.

## BODYFLEX II



Cooler and more  
Comfortable.

ITEM # 45-729	MEASURE CENTER OF KNEE
45-729-XS	13" TO 14"
45-729-S	14" TO 15"
45-729-M	15" TO 16"
45-729-L	16" TO 17"
45-729-XL	17" TO 19"
45-729-XXL	19" TO 21"