

CINDY WRIST BRACE

FEATURES

- Flexible palmar and dorsal stays help maintain wrist in a neutral position.
- Allow full dexterity of the fingers and thumb, yet limits wrist flexion.
- Constructed of durable, breathable, padded material.
- Loop and lock closure for easy adjustment.
- 6" in length.



BAG WITH INSERT

ITEM # LEFT HAND	ITEM # RIGHT HAND
BC5541-LXS	BC5541-RXS
BC5541-LS	BC5541-RS
BC5541-LM	BC5541-RM
BC5541-LL	BC5541-RL
BC5541-LXL	BC5541-RXL

WRIST CIRCUMFERENCE

X-SMALL (UP TO 5 1/2")
 SMALL (5 1/2" - 6 1/2")
 MEDIUM (6 1/2" - 7 1/2")
 LARGE (7 1/2" - 8 1/2")
 X-LARGE (8 1/2" - 9 1/2")

RETAIL PACKAGING

ITEM # LEFT HAND	ITEM # RIGHT HAND
PC5541-LXS	PC5541-RXS
PC5541-LS	PC5541-RS
PC5541-LM	PC5541-RM
PC5541-LL	PC5541-RL
PC5541-LXL	PC5541-RXL

WRIST CIRCUMFERENCE

X-SMALL (UP TO 5 1/2")
 SMALL (5 1/2" - 6 1/2")
 MEDIUM (6 1/2" - 7 1/2")
 LARGE (7 1/2" - 8 1/2")
 X-LARGE (8 1/2" - 9 1/2")