

# F8 ANKLE BRACE

## FEATURES

- Designed to support and protect ankles.
- Figure 8 configuration eliminates taping.
- Thin yet strong ballistic type nylon is less bulky for a comfortable fit.
- Figure 8 straps are non-stretch lightweight nylon.
- Elastic panel at the Achilles provides superior comfort and fit.
- Knit tongue is cool and provides for easy application.



ITEM # BAG WITH INSERT	ITEM # RETAIL PACKAGED	SHOW SIZE
BCF8-XS	PCF8-XS	MALE (3 - 5) FEMALE (4-6)
BCF8-S	PCF8-S	MALE (6-7) FEMALE (7-9)
BCF8-M	PCF8-M	MALE (8-10) FEMALE (10-11)
BCF8-L	PCF8-L	MALE (11-12) FEMALE (12-13)
BCF8-XL	PCF8-XL	MALE (13-14) FEMALE (14+)
BCF8-XXL	PCF8-XXL	MALE (15+)