

HUMERAL BRACE SHOULDER

FEATURES

- Treatment of diaphyseal fractures.
- Deltoid cap provides greater alignment control & prevents distal slippage.
- Complete with two double thickness cotton stockinettes.
- Foam lined clamshell providing superior compression to the fracture sight.



ITEM #	BICEP CIRCUMFERENCE	BICEP LENGTH
HFODE-S	8" TO 11"	11"
HFODE-M	11" - 14"	12.5"
HFODE-L	12" - 15"	13.5"
HFODE-XL	15" - 16"	14"