

# K15

## FEATURES

- 1/8" nylon 2-sides neoprene.
- Patellar opening for positioning & comfort.
- Patellar pull-panel w/ adjustable straps.
- 1/4" tubular, enclosed j-pad buttress.
- Medial & Lateral spiral stays.



ITEM # LEFT	ITEM # RIGHT
K15-LXS	K15-RXS
K15-LS	K15-RS
K15-LM	K15-RM
K15-LL	K15-RL
K15-LXL	K15-RXL
K15-LXXL	K15-RXXL

MEASURE CENTER OF KNEE
12" TO 13"
13" TO 14"
14" TO 15"
15" TO 17"
16" TO 18"
18" TO 20"