

# K64

## FEATURES

- Indications: Arthritis, Post-Operative management, Anterior & Posterior cruciate instability, Medial and Lateral Collateral instability.
- 3/16 nylon 2-sides neoprene.
- Patellar opening for positioning and comfort.
- Upper & Lower circumferential straps w/ finger loops & pull tabs.
- Adjustable elastic segment to prevent tourniquet effect in popliteal opening.
- Large, comfortable popliteal opening.



ITEM #	MEASURE CENTER OF KNEE
K64-XS	12" TO 13"
K64-S	13" TO 14"
K64-M	14" TO 15"
K64-L	15" TO 17"
K64-XL	16" TO 18"
K64-XXL	18" TO 20"