


KNEE COMPRESSION WITH STAYS

FEATURES

- Easy slip on, soothing compression, open or closed patella, latex free
- Unique 4-way stretch
- Breathable material is excellent alternative to neoprene
- Seamless construction to prevent irritation
- Knee strains, Chondromalacia, Swelling & Edema, Post Operative
- Recurrent knee pain, tendonitis, Arthritis, Protection
- Spiral, stainless steel stays for a higher level of support



ITEM # TAN CLOSED	ITEM # TAN OPEN	4.5" BELOW KNEE	5.5" ABOVE KNEE
300T-XS	300POT-XS	9.75" TO 11"	14.5" TO 16"
300T-S	300POT-S	11" TO 12"	16" TO 17"
300T-M	300POT-M	11.5" TO 13.25"	17" TO 18"
300T-L	300POT-L	13.25" TO 14.25"	18" TO 19.25"
300T-XL	300POT-XL	14.25" TO 15.5"	19.25" TO 20.5"
300T-XXL	300POT-XXL	15.5" TO 16.5"	20.5" TO 22"
300T-XXXL	300POT-XXXL	16.5" TO 18"	22" TO 23"

ITEM # BLACK CLOSED	ITEM # BLACK OPEN	4.5" BELOW KNEE	5.5" ABOVE KNEE
300B-XS	300POB-XS	9.75" TO 11"	14.5" TO 16"
300B-S	300POB-S	11" TO 12"	16" TO 17"
300B-M	300POB-M	11.5" TO 13.25"	17" TO 18"
300B-L	300POB-L	13.25" TO 14.25"	18" TO 19.25"
300B-XL	300POB-XL	14.25" TO 15.5"	19.25" TO 20.5"
300B-XXL	300POB-XXL	15.5" TO 16.5"	20.5" TO 22"
300B-XXXL	300POB-XXXL	16.5" TO 18"	22" TO 23"