

Air Walker

Designed for orthopaedic foot and ankle injuries



Features*



Dual air chambers positioned medial and lateral at the ankle

Three-piece design includes posterior and anterior shells for true "cast-like" immobilization

Anti-migration strap to prevent forefoot from sliding

Deluxe liner with detachable toe cover

* Patent Pending

Ergo-dynamic, lightweight design in modern colour | Wide, low-profile rocker sole to facilitate natural gait and provide increased stability | Wider foot base for improved comfort and support | Economically priced | Latex-free

Order Information

Reorder No.	Size	Male Shoe	Female Shoe
BC0982	Small	4-7	5-8
BC0983	Medium	7 ½ -11	8 ½ -12
BC0984	Large/X-Large	11 ½ +	12 ½ +

Landmark Medical Systems

#19-100 West Beaver Creek | Richmond Hill | Ontario | Canada | L4B 1H4

T 905.881.2324 | F 905.881.2327