

PREMIER WITH THUMB

FEATURES

- Indicated when immobilization of the MP & CMC joint is desired.
- Maximum contouring at the palmar for full finger function.
- Removable and adjustable palmar stay for maximum support.
- Removable rigid thumb stay is adjustable for custom fit.
- Constructed of durable, breathable, padded material.
- Loop and lock closure for easy adjustment.
- 8" in length



BAG WITH INSERT

ITEM # LEFT HAND	ITEM # RIGHT HAND
BC4551-LXS	BC4551-RXS
BC4551-LS	BC4551-RS
BC4551-LM	BC4551-RM
BC4551-LL	BC4551-RL
BC4551-LXL	BC4551-RXL

WRIST CIRCUMFERENCE
X-SMALL (UP TO - 6 1/2")
SMALL (6 1/2" - 7 1/4")
MEDIUM (7 1/4" - 8 1/4")
LARGE (8 1/4" - 9")
X-LARGE (9" - 10")

RETAIL PACKAGING

ITEM # LEFT HAND	ITEM # RIGHT HAND
PC4551-LXS	PC4551-RXS
PC4551-LS	PC4551-RS
PC4551-LM	PC4551-RM
PC4551-LL	PC4551-RL
PC4551-LXL	PC4551-RXL

WRIST CIRCUMFERENCE
X-SMALL (UP TO - 6 1/2")
SMALL (6 1/2" - 7 1/4")
MEDIUM (7 1/4" - 8 1/4")
LARGE (8 1/4" - 9")
X-LARGE (9" - 10")