

# VACO ANKLE BRACE

## FEATURES

- Provides firm support for upper ankle joint against pro and supination by means of a construction consisting of a medial and lateral shell
- Shells are interconnected by a clamp system for easy adjustment
- Shells can be moved ventrally to ensure firm stability and promote comfort
- A vacuum cushion allows pressure-free adaptation to address changing swelling conditions



ITEM # VACO ANKLE	HEIGHT IN INCHES
43-VANK1	HEIGHT GREATER THEN 5'3"
43-VANK2	HEIGHT OF 5'3" OR SHORTER